





UPLC/UHPLC Cortecs kolóny od WATERS

Rozšírite možnosti Vášho analytického laboratória s najnovším prírastkom chromatografických kolón od spoločnosti Waters. CORTECS® 1.6 µm a 2.7 µm kolóny s pevným jadrom ponúkajú novú úroveň separačnej účinnosti a rýchlejšie analýzy v porovnaní s plne poréznymi časticami. V ponuke nájdete UPLC ako aj HPLC formáty v prevedení C18, C18+ a HILIC. V ponuke sú aj predkolónky s rovnakými druhmi sorbentov. Viac informácií nájdete na oficiálnych stránkach výrobcu.

Schéma znázorňujúca dostupné formáty pre Cortecs kolóny.



CORTECS[™]
COLUMNS

	C ₁₈ +	C ₁₈	HILIC
Chemistry			
Intended Use	General purpose, high-efficiency, reversed-phase column. A positively charged surface delivers excellent peak shape for basic compounds at low pH.	General purpose, high-efficiency, reversed-phase column. Balanced retention of acids, bases and neutrals at low and mid-range pH.	High-efficiency column designed for retention of extremely polar analytes. Offers orthogonal selectivity vs. C ₁₈ columns.

 **LAMBDA LIFE**